

# Ellie Recliner

## ELIXXXXXX

kwalu.

### Unite Style



**1**  
23" SW



**2**  
Trendelenburg

### Arm Style



**1**  
Open Loop Arm



**2**  
Upholstered Arm Panels


### Legs Style



**A**  
Locking Casters



**B**  
Central Locking System



**C**  
Protea™ Legs

### \*Transfer Arm



**0**  
No Transfer




**1**  
Left Transfer




**2**  
Right Transfer

### Pillow




**0**  
No Pillow




**1**  
With Pillow

### Lever Position



**R**  
Right Arm



**L**  
Left Arm

### Available Options



**\*\*OPT#48**  
IV Pole Holder



**\*\*OPT#189**  
Foley Hook



**\*\*\*OPT#201**  
Heat Massage

### \*\*\*\*Table Options



**OPT#57EF**  
Folding Side Table



**OPT#57ESS**  
Tablet Table Solid Surface



**OPT#57ESP**  
Tablet Table Protea Fusion

### Notes

\*Transfer Arm and Lever cannot be placed on the same arm.

\*\*Must be placed on the opposite side of lever.

\*\*\*Available only when Upholstered Arm Panels is selected. Heat Massage controller will always be located on same side as lever.

\*\*\*\*The Tablet and Side Folding Tables must be positioned on the opposite side of the lever and cannot be used with the Transfer Arm option.