

### Kwalu Features

# Product Dimensions and Specs

 Overall
 29W 31D 35H

 Seat
 23W 18D 20H

Arm 26H Wt. 49 lbs Wt. Capacity 350 lbs

# **Product Configurations**

Pisa Collection	PISAXXXXX
Number of Seats	1 Single Love Seat
Pillow	With Pillow No Pillow
Back Detail	1 No Detail With Channel Stitching
Leg Style	A Tapered
Seat Construction	LS Lose Seat Cushion (Recommended for Senior Living application)  Fixed Seat Cushion (Recommended for Healthcare application)

Note

<sup>\*</sup> Overall Depth will be shorter if this option is selected



### Product Statement of Line Images



Pisa Lounge Chair, Fixed Seat Pisa Lounge Chair, Loose Model: PISA1A1AFS 29W 31D 35H 26AH



Seat Cushion Model: PISA1A1ALS 29W 31D 35H 26AH



Pisa Lounge Chair, with Pillows, with Channel Stitching Fixed Seat Model: PISA1A2AFS 29W 31D 35H 26AH



Pisa Lounge Chair, Loose Seat Cushion, with Pillows, with Channel Stitching Model: PISA1A2ALS 29W 31D 35H 26AH



Pisa Lounge Chair, without Pillows, without Channel Stitching Fixed Seat Model: PISA1B1AFS 29W 26D 30H 26AH



Pisa Lounge Chair, Loose Seat Cushion, without Pillows, without Channel Stitching Model: PISA1B1ALS 29W 26D 30H 26AH




Pisa Lounge Chair, without Pillows, with Channel Stitching Fixed Seat Model: PISA1B2AFS 29W 26D 30H 26AH




Pisa Lounge Chair, Loose Seat Cushion, without Pillows, without Channel Stitching Model: PISA1B2ALS 29W 26D 30H 26AH







Pisa Love Seat, with Pillows, without Channel Stitching Fixed Seat Model: PISA2A1AFS 60W 31D 35H 26AH



Pisa Love Seat, Loose Seat Cushion Model: PISA2A1ALS 60W 31D 35H 26AH



Pisa Love Seat, with Pillows, with Channel Stitching Fixed Seat Model: PISA2A2AFS 60W 31D 35H 26AH




Pisa Love Seat, Loose Seat Cushion, with Pillows, with Channel Stitching Model: PISA2A2ALS 60W 31D 35H 26AH




Pisa Love Seat, without Pillows, without Channel Stitching Fixed Seat Model: PISA2B1AFS 60W 26D 30H 26AH



Pisa Love Seat, Loose Seat Cushion, without Pillows, without Channel Stitching Model: PISA2B1ALS 60W 26D 30H 26AH



Pisa Love Seat, without Pillows, with Channel Stitching Cushion, without Pillows, Fixed Seat Model: PISA2B2AFS 60W 26D 30H 26AH



Pisa Love Seat, Loose Seat without Channel Stitching Model: PISA2B2ALS 60W 26D 30H 26AH



#### **Product Features**



Sculpted Leg Detail

Sculpted Leg Detail enhances the look of the product making it more elegant.



Steel Reinforced

Our patented steelreinforced joint system uses over 20 pieces of steel in most Kwalu chairs. Stop ongoing and costly furniture replacement.

<u>Learn More ></u>



Attached Back Cushion

Back cushion is attached with D-rings which allows the caregiver to remove the cushion when needed.



Removable Seat Deck

Removable seat deck prevents buildup of debris and provides for easy maintenance.



Number of Seats

Choose from 2 seat styles - lounge chair and love seat.

<u>Learn more ></u>



Optional Seat

Choose the best seat cushion for your application, a non-removable seat cushion for healthcare, or a loose seat cushion for senior living.

#### **Product Foam**

Foam is compliant with CAL 117 and AB 2998

#### **Product Glides**

Clear Plastic Butyl Screw in Glide



### **Product Finishes**

#### Base Finish





## Product Environmental Certifications



